

IEC 60601-1 Medical electrical equipment

Part 1: General requirements for basic safety and essential performance

 Report Reference No.......
 MDT2013-0016

 Date of issue
 2014-11-27

 Total number of pages
 194 pages

Tested by (name + signature): KI IL NAM

Approved by (+ signature)...... Hong-Soon Choi

Testing Laboratory...... KTC (Korea Testing Certification)

Address 22 Heungan-daero 27 beon-gil, Gunpo-si, Gyeonggi-Do, 435-862

KOREA REPUBLIC OF

Applicant's name...... R&L Co., Ltd.

Address (Soha-dong, 11 Floor, KwangMyung Ace High-End Tower,

Sohadong),108, Haan-ro, Gwangmyeong-si, Gyeonggi-do, Korea

Test specification:

Standard IEC 60601-1: 2005 + CORR. 1 (2006) + CORR. 2 (2007)

Test procedure...... Test Report

Non-standard test N/A

method....:

Test Report Form No...... IEC60601 1G

Test Report Form Originator......: Underwriters Laboratories Inc.

Master TRF Dated 2011-11

Copyright © 2011 IEC System for Conformity Testing and Certification of Electrical Equipment (IECEE), Geneva, Switzerland. All rights reserved.

This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.

If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo shall be removed

This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.

Test item description BIO MAT

Trade Mark:

Manufacturer...... R&L Co., Ltd.

Model/Type reference...... BIO MAT PROFESSIONAL S

Ratings 230 V~, 50 Hz, 200 W

List of Attachments (including a total number of pages in each attachment):	
1. Attachment #1 : Photographs (Page 190 ~ page 194 in this report)	
Summary of testing	
Tests performed (name of test and test clause):	Testing location:
Refer each Clause section	KTC (Korea Testing
All applicable test performed and all applicable clauses are evaluated with a exception of Usability(Clause 7.1.1 & 12.2)	Certification)
	22 Heungan-daero 27 beon-gil, Gunpo-si, Gyeonggi-Do, 435-862 KOREA REPUBLIC OF

Page 2 of 194